



ME 5 & ME 10

Efficient, lightweight, ME5 and ME10 are indicated for measuring functional oxygen saturation of arterial hemoglobin (SpO2) and pulse rate of adult and pediatric patients (>25 kg).

Technical Specifications





ME 5

Apt to spot checking of adult and pediatric patients (>25kg) on finger, ME5 displays real-time SpO2 measurements, pulse, and the quality of the signal.

SpO2 range: 35 - 99%,Pulse rate: 30 - 235 BPM

- Incl. bag, neck string and 2 x AAA batteries

- Dimensions: L 58 x W 32 x H 34 mm

- Weight: 50 g (with batteries)

ME 10

ME10 pulse oximeter has a color LCD screen, 6 display options, plethysmographic wave display, SpO2 and heart rate, visual alarm of the battery level and auto-off function after 8 seconds of non-use.

- Switchable display in 6 different positions
- Color and numeric display of SpO2, Pulse and Plethysmographic wave
- SpO2 range: 45 100%,
- Pulse rate: 20 300 BPM
- Incl. bag, neck string, 2x AAA batteries
- Dimensions: L 60 x W 34 x H 30 mm
- Weight: 50g (with batteries)



Im Erlengrund 20, 46149 Oberhausen - Germany

t: +49 (0) 208 377 890-0 info@medical-econet.com f: +49 (0) 208 377 890-55 www.medical-econet.com facebook.com/medicaleconetgermany twitter.com/medicaleconet

٦i	cti	rik	NI I	tο	r